



Medical Catalog for 2023-2024





Aeternus Medical operates with a single-point-of-contact model, managing the triad of devices, consumables, and services end-to-end. We begin with a site survey and needs analysis, plan the project with 2D layouts and infrastructure guidance, and then ship and install the solutions you select.

Device portfolio: BioSTI & XGATE implant systems; intraoral scanners; Formlabs Dental 3D printers; dental units and dentist/assistant stools; compressors & central vacuum, the sterilization chain (autoclave–ultrasonic–pouch), RVG/Panoramic/CBCT imaging (optional), and CAD/CAM infrastructure.

Consumables & supplies: clinical disposables (masks–gloves–drapes, suction tips, irrigation sets, disinfectants); restorative materials (bonding–composites–cements); endodontic consumables and irrigants; impression materials–plaster–wax; laboratory supplies (zirconia/PMMA blocks, burs, discs, polishing stones/rubbers, sandblasting–polishing media, furnace accessories) and maintenance chemicals.

Process & service: After delivery and installation, we perform IQ/OQ tests, provide user training, and hand over as-built documentation; we ensure continuity with warranty, spare parts, and periodic maintenance plans. For consumables, we offer steady stock, fast delivery, and contracted pricing advantages.

Aeternus drives the digital transformation of clinics and laboratories with a transparent timeline, measurable performance, and strong after-sales support. Our goal is to accelerate clinical workflows, make costs predictable, and achieve repeatable clinical success for every patient.

Index

01-06	LED Low speed handpiece
07-12	High speed handpiece
13-20	Ultrasonic scaler
21-26	Air scaler
27-30	Teeth whitening machine
31-32	Teeth whitening accessories
33-38	Root Apex locator
39-42	Endo motor
43-45	Cutter gutter
46	Endodontic equipment
47-48	Obturation system
49-54	Dental curing light
55-56	Dental Oral Microscope
57-66	LED Dental Implanting Lamp
67-68	DENTAL Surgical loupes
69-72	Dental unit
73-76	Portable dental unit
77-78	Compressor
79-84	Dental cabinet
85-86	Dentist chair
87-90	Micromotor
91-93	Intra camera
95-102	X ray machine and sensor
103-104	Steam Sterilizer
105-106	Ultrasonic Cleaner
107-108	Dental lab equipment
109-110	Dental tool
111-114	Dental disposables

20:1 Implant contra angle



Order code.:AM-SG20

Chuck type push button
Bur applicable $\Phi 2.35\text{mm}$
Rotation speed $\leq 0 \sim 2.500\text{rpm}$
Max torque 80N.M
Noise $\leq 65\text{db}$
Ratio:20:1

1:5 Increasing Contra Angle



Order code.:AM-Z95L

Chuck type push button
Bur applicable $\Phi 1.60\text{mm}$
Rotation speed $\leq 200.000\text{ rpm}$
Spray 4 water spray
Air pressure 0.22-0.25MPa
Noise $\leq 68\text{db}$
Ratio:1:5

20:1 Fiber optic contra angle



Order code.:AM-SG20F

Chuck type push button
Bur applicable $\Phi 2.35\text{mm}$
Rotation speed $\leq 0 \sim 2.500\text{rpm}$
Max torque 80N.M
Noise $\leq 65\text{db}$
Ratio:20:1

1:4.2 Increasing Contra Angle



Order code.:AM-Z45L

Low speed handpiece 1:4.2 speed
Chuck type push button
Bur applicable $\Phi 1.60\text{mm}$
Rotation speed $\leq 168.000\text{ RPM}$
Spray Inner channel
Noise $\leq 68\text{db}$
Ratio:1:4.2

20:1 LED contra angle



Order code.:AM-SG20L

Chuck type push button
Bur applicable $\Phi 2.35\text{mm}$
Rotation speed $\leq 0 \sim 2.500\text{rpm}$
Max torque 80N.M
Noise $\leq 65\text{db}$
Ratio:20:1

Low speed handpiece set



Order code:AM-K41



Rotation speed:1:1
Chuck type:Push button
Bur applicable: $\Phi 2.35\text{mm}$
Rotation speed:40,000 rpm
Air pressure:0.22-0.25MPa
Noise $\leq 65\text{db}$



Order code:AM-K42

Low-speed Handpiece Set



Order code:AM-N44



FPB-B



AM-43S



AM-N43M



Order code:AM-N43



Order code:AM-N41



Order code:AM-N44



AM-FX23



AM-FX65



AM-FX205



Order code : AM-BB



With double 2.78 mm
ceramic bearing,
rotation controlling



Order code : AM-FPB

Low speed handpiece



Order code:AM-N42-C



Order code:AM-N42-D



Order code:AM-N42-D



Order code:AM-RA41

64:1 Endodontic contra angle



Order code : AM-ER64

Chuck type push button

Bur applicable $\Phi 2.35\text{mm}$

Rotation speed $\leq 625\text{ rpm}$

Spray external

Rotating way circle

10:1 Endodontic contra angle



Order code : AM-ER10

Chuck type push button

Bur applicable $\Phi 2.35\text{mm}$

Rotation speed $\leq 4000\text{ rpm}$

Spray external

Rotating way 30-30 degree rotating

Surgical straight Handpiece




Order code : AM-SH01



Saw Handpieces



Order code : AM-SH02



High Speed Handpiece



Order code:AM-Max-L2

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	3 water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-Max

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	3 water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-Max3L

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	4 water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-T3

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	4 water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-98L1

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	Single water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-5LED

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	5 water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-95L

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	Single water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-SMD

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	Four water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-636L

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	3 water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-T3L

Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	Single water spray
Air pressure	0.22-0.25MPa
Noise	≤68db

High Speed Handpiece



Order code:AM-636L

Chuck type Push button
Bur applicable $\Phi 1.58-1.6\text{mm}$
Rotation speed $\geq 360.000\text{ rpm}$
Spray 3 water spray
Air pressure $0.22-0.25\text{MPa}$
Noise $\leq 68\text{db}$



Order code:AM-45L

Chuck type Push button
Bur applicable $\Phi 1.58-1.6\text{mm}$
Rotation speed $\geq 360.000\text{ rpm}$
Spray 4 water spray
Air pressure $0.22-0.25\text{MPa}$
Noise $\leq 68\text{db}$



Order code:AM-45

Chuck type Push button
Bur applicable $\Phi 1.58-1.6\text{mm}$
Rotation speed $\geq 360.000\text{ rpm}$
Spray 4 water spray
Air pressure $0.22-0.25\text{MPa}$
Noise $\leq 68\text{db}$



Order code:AM-98L1

Chuck type Push button
Bur applicable $\Phi 1.58-1.6\text{mm}$
Rotation speed $\geq 360.000\text{ rpm}$
Spray Single water spray
Air pressure $0.22-0.25\text{MPa}$
Noise $\leq 68\text{db}$



Order code:AM-MN02

Chuck type Push button
Bur applicable $\Phi 1.58-1.6\text{mm}$
Rotation speed $\geq 360.000\text{ rpm}$
Spray Single water spray
Air pressure $0.22-0.25\text{MPa}$
Noise $\leq 68\text{db}$



Order code:AM-SMD

Chuck type Push button
Bur applicable $\Phi 1.58-1.6\text{mm}$
Rotation speed $\geq 360.000\text{ rpm}$
Spray Four water spray
Air pressure $0.22-0.25\text{MPa}$
Noise $\leq 68\text{db}$

K series Quick coupling handpiece



Order code:AM-8000B



Order code:AM-8000M



Order code:AM-9000L



Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	4 water spray
Air pressure	0.22-0.25MPa
Noise	≤68db



Order code:AM-X600L



Order code:AM-X600LN



Order code:AM-X600LQ



Order code:AM-X600LN-Q



Chuck type	Push button
Bur applicable	Φ1.58-1.6mm
Rotation speed	≥360.000 rpm
Spray	4 water spray
Air pressure	0.22-0.25MPa
Noise	≤68db

Ultrasonic scaler

DTE D7 ultrasonic scaler



Parameter:

Function:Scaling,Perio,Endo
Handpiece:LED handpiece
Power supply Input: 220-240V~ 50Hz/60Hz
Main unit input: 24V~ 1.3A
Output primary tip vibration excursion: $\leq 100\mu\text{m}$
Output tip vibration frequency: $28\text{kHz}\pm 3\text{kHz}$
Output power: 3W to 20W
Weight of main unit: 1.17kg

Order code:AM-001

功能：洁治、牙周、根管荡洗
工作手柄：发光手柄
工作电压：220-240V~ 50Hz/60Hz
功率：3W-20W
震动频率：28kHz $\pm$ 3kHz

DTE D600 Ultrasonic scaler



Parameter:

Function:Scaling,Perio,Endo
Handpiece:LED handpiece
Main unit input: 220-240V~ 50Hz/60Hz 150mA
Output primary tip vibration excursion: $\leq 90\mu\text{m}$
Output half-excursion force: $< 2\text{N}$
Output tip vibration frequency: $28\text{kHz}\pm 3\text{kHz}$
Output power: 3W to 20W
Main unit fuse: T0.5AL 250V
Water pressure: 0.01MPa to 0.5MP

Order code:AM-003

功能：洁治、牙周、根管荡洗
工作手柄：发光手柄
工作电压：220-240V~ 50Hz/60Hz
功率：3W-20W
震动频率：28kHz $\pm$ 3kHz

UDS-E Ultrasonic scaler



Parameter:

Function:Scaling,Perio,Endo
Handpiece:LED handpiece
Power supply Input: 220-240V~ 50Hz/60Hz
Main unit input: 24V~ 1.3A
Output primary tip vibration excursion: $\leq 100\mu\text{m}$
Output tip vibration frequency: $28\text{kHz}\pm 3\text{kHz}$
Output power: 3W to 20W
Weight of main unit: 1.17kg

Order code:AM-002

功能：洁治、牙周、根管荡洗
工作手柄：发光手柄
工作电压：220-240V~ 50Hz/60Hz
功率：3W-20W
震动频率：28kHz $\pm$ 3kHz

U600 ultrasonic scaler



Parameter:

Function:Scaling,Perio,Endo
Handpiece:LED handpiece
Main unit input: 220-240V~ 50Hz/60Hz 150mA
Output primary tip vibration excursion: $\leq 90\mu\text{m}$
Output half-excursion force: $< 2\text{N}$
Output tip vibration frequency: $28\text{kHz}\pm 3\text{kHz}$
Output power: 3W to 20W
Main unit fuse: T0.5AL 250V
Water pressure: 0.01MPa to 0.5MP

Order code:AM-004

功能：洁治、牙周、根管荡洗
工作手柄：发光手柄
工作电压：220-240V~ 50Hz/60Hz
功率：3W-20W
震动频率：28kHz $\pm$ 3kHz

Ultrasonic scaler

VRN A8 Ultrasonic scaler



Function : Scaling,Perio,Endo
Handpiece : Detachable handpiece
Foot switch battery : 1.5V*2
Receiving sensitivity : -114dB
Water pressure : 0.01Mpa-0.5Mpa
Input : AC100~240V 47~63Hz
Power : 30VA-48VA
PC : 30V 1.5A

Vibration frequency of the tip : <25~31KHz

Order code:AM-005

功能：洁治、牙周、根管荡洗
工作手柄：发光手柄
工作电压：220~240V~ 50Hz/60Hz
功率：3W~20W
震动频率：25kHz±3kHz

Vibration excursion of the tip : <100μM
Ultrasonic output power : 3W~20W
Weight of main unit : 1.1KG
G.W : 2.50KG
Dimensions : 330 x 155 x 155 mm
Packing size : 340 x 240 x 180mm

VRN Q5 Ultrasonic scaler



Function : Scaling,Perio,Endo
Handpiece : Detachable handpiece
Foot switch battery : 1.5V*2
Receiving sensitivity : -114dB
Water pressure : 0.01Mpa-0.5Mpa
Input : AC100~240V 47~63Hz
Power : 30VA-48VA
PC : 30V 1.5A
Vibration frequency of the tip : <25~31KHz
Vibration excursion of the tip : <100μM
Ultrasonic output power : 3W~20W

Order code:AM-007

功能：洁治、牙周、根管荡洗
工作手柄：发光手柄
工作电压：220~240V~ 50Hz/60Hz
功率：3W~20W
震动频率：25kHz±3kHz

VRN-B ultrasonic scaler



Host input DC30V 1.5A
Working handle Siamese handle
Tip vibration 25~31KHz Primary
Function cleansing,
Tip main vibration ≤100um
Control method 3W~20W
Working water 0.01MPa~0.05MPa
Power input AC100~240V 47~63Hz

Working voltage 30VA~48VA

Order code:AM-006

功能：洁治、牙周、根管荡洗
工作手柄：不发光手柄
工作电压：220~240V~ 50Hz/60Hz
功率：3W~20W
震动频率：25kHz±3kHz

VRN Q6 Ultrasonic periodontal Therapy System



Main Unit Input : DC30V 1.3A
Water Input pressure : 0.5Mpa~0.6Mpa
Foot Switch Batteries : 1.5V*2
Input voltage : 30VA~48VA
Ultrasonic output power : 3W~20W
Weight of main unit : 1.5kg
Gross Weight : 6kg
Main Unit Dimensions : 310mm x 370mm x170mm

Order code:AM-008

主机输入：DC30V 1.3A
运行模式：连续运行
输入功率：30VA~48VA
超声输出功率Ultrasonic output power: 3W~20W
包装尺寸Packing Size: 470×400×235mm
熔断器：T1AL 250V
喷砂粉规格：龈上喷砂粉为100目~200目
龈下喷砂粉应>200目，微咸或柠檬味

Ultrasonic scaler

B5 ultrasonic scaler



Order code:AM-009

手柄型号: H1
震动频率: 25~31kHz
适配电压: 100~240v
包装尺寸: 25.5*21.5*19.5cm
包装重量: 1.75kg

Function:Scaling,Perio
Handpiece material: Imported plastic
Handpiece color: Light blue
Autoclavable: Not allowed
Certificate: CE FDA ISO
Usage: By dentist or veterinary
Voltage: 100-240V
Instrument classification: Class II
Tip vibration frequency: 28kHz±3kHz
Handpiece item:H1

P5 ultrasonic scaler



Order code:AM-010

手柄型号: H3
震动频率: 25~31kHz
适配电压: 100~240v
包装尺寸: 25.5*21.5*19.5cm
包装重量: 1.75kg

Function:Scaling,Perio
Main unit input:30V 1.2A
Handpiece color: Light blue
Output power:3W-20W
Output tip vibration frequency:28kHz±3kHz
Water pressure:0.1Mpa~0.5MPa (1bar~5bar)
Voltage: 100-240V
Weight of main unit: 0.55kg

Max Piezo 7+ ultrasonic scaler



Order code:AM-011

功能: 洁治、牙周、根管荡洗
工作手柄: 发光手柄
工作电压: 220~240V~ 50Hz/60Hz
震动频率: 28kHz±3kHz
功率: 38W
尖端主震荡偏移: < 90 μm

Equipment function: G、P、E
(With root canal washing function)
Frequency: 28kHz±3kHz
Output power: 38VA
Tip main vibration deviation: < 90μm
Tip output characteristics:semi-offset force0.5 N~2N
Power Input:220V 50Hz/110V 50Hz
Water pressure: 0.01MPa~0.5MPa
Water supply method:Self water supply
Tip output power:3W~20W
Host insurance: T0.5AL 250VAC

DS7+ Ultrasonic scaler



Order code:AM-012

功能: 洁治、牙周、根管荡洗
工作手柄: 发光手柄
工作电压: 220~240V~ 50Hz/60Hz
震动频率: 28kHz±3kHz
功率: 38W
尖端主震荡偏移: < 90 μm

Equipment function: G、P、E
(With root canal washing function)
Frequency: 28kHz±3kHz
Output power: 38VA
Tip main vibration deviation: < 90μm
Tip output characteristics:semi-offset force0.5 N~2N
Power Input:220V 50Hz/110V 50Hz
Water pressure: 0.01MPa~0.5MPa
Water supply method:Self water supply
Tip output power:3W~20W
Host insurance: T0.5AL 250VAC

Ultrasonic scaler

AP-A Dental Air Polisher



Order code:AM-013

工作电压: 100-240V~ 50/60Hz
主机输入: 30V 1.3A
主机保险: T1.6AL 250V
进气压力: 5.5bar~7.5bar (0.55MPa~0.75MPa)
主机重量: 2.0Kg
主机尺寸: 270mm×170mm×90mm

Power adapter input: 100-240V~ 50/60Hz
Main unit input: 30V 1.3A
Main unit insurance: T1.6AL 250V
Air inlet pressure: 5.5bar ~ 7.5bar(0.55MPa ~ 0.75MPa)
Main unit weight: 2.0kg
Main unit size: 270mm×170mm×90mm

Refine PTX Dental Scaler and Air Polisher



Order code:AM-015

工作电压: 100-240V~ 50/60Hz
主机输入: 30V 1.3A
主机保险: T1.6AL 250V
进气压力: 5.5bar~7.5bar (0.55MPa~0.75MPa)
主机重量: 4.0Kg
主机尺寸: 320*330*120mm

Power input:100-240V~ 50/60Hz
Main unit input: 30V 1.3A
Control mode:Touch Control
Water supply Mode:Water supply
Power Adaptor Input:100-240V~ 50/60Hz
Main unit fuse:T1.6AL 250V
Water inlet pressure:0.1bar-5bar
Air inlet pressure:5.5bar-7.5bar
Main unit Size:320*330*120mm
Main unit weight:4.0kgs

AP-A Dental Scaler and Air Polisher



Order code:AM-014

工作电压: 100-240V~ 50/60Hz
主机输入: 30V 1.3A
输出的尖端主振动偏移(最大值): 100 μm, 偏差: ±50%
输出的尖端振动频率: 30±5kHz
尖端输出功率: 3W~20W
主机保险: T1.6AL 250V

Power adapter input: 100-240V~ 50/60Hz 1.1A
Main unit input: 30V 1.3A
Main vibration offset of output tip (maximum):100μm, deviation: ± 50%
Output tip vibration frequency: 30±5kHz
Output half offset force (maximum):10N deviation: ± 50%
Output power of tip: 3W ~ 20W
Main unit insurance: T1.6AL 250V
Air inlet pressure: 5.5bar ~ 7.5bar (0.55MPa ~ 0.75MPa)
Main unit weight:2.0kg
Main unit size: 270mm×170mm×90mm

2 in 1 Complex Scaler



Order code:AM-016

输入功率 8W
电压 AC100V~240V 50~60Hz
尖端主振动偏移 0.05-0.15 μm,
工作尖输出频率
偏移功率 ≥0.5N
环境温度 5℃-40℃
相对湿度 ≤80%
气压 86kpa-106kPa

Input Power 8W
Voltage AC100V ~ 240V 50 ~ 60Hz
Transformer: DC24V
Main libration excursion 0.05-0.15μm ,
Tips output libration freque 25KHz ~ 35KHz的
Operation type Continue operation
Environment Temperature 5℃-40℃
Relative humidity ≤80%
Atmospheric pressure 86kpa-106kPa

Ultrasonic scaler

UDS-N2 Built-in Ultrasonic Scaler



Order code:AM-017

UDS-N3 LED Built-in Ultrasonic Scaler



Order code:AM-018

Scaler handpiece



Order code:AM-019

Scaler working tips



Order code:AM-020

Air Scaler



Order code:AM-021

Air Scaler



Order code:AM-022

Air Scaler



Order code:AM-023

Air Scaler



Order code:AM-024

Ultrasonic scaler

Gingival Flusher



Order code:AM-025

Air polisher



Order code:AM-029

Quick coupling Air polisher



Order code:AM-033

VRN Scaler Automatic Water Feeder



Order code:AM-026

Air polisher



Order code:AM-030

Air polisher



Order code:AM-034

Scaler Automatic Water Feeder



Order code:AM-027

Luxury Air polisher



Order code:AM-031

Air polisher



Order code:AM-035

Scaler Automatic Water Feeder



Order code:AM-028

Coxo Air polisher



Order code:AM-032

Polish Sand powder



(Lemon/Mint/Strawberry)

Order code:AM-036



(Lemon/Mint/Strawberry)

Order code:AM-037

Ultrasonic scaler

Alumina Air Polisher



Order code:AM-038

ART-M1 Magnetic Stretch Type
Ultrasonic Scaler



Order code:AM-042

Dental dust box



Order code:AM-046

Alumina Air Polisher



Order code:AM-039

ATR-M3 II Magnetostrictive painless
periodontal treatment Instrument



Order code:AM-043

Dental Dust Collector



Order code:AM-047

Alumina Air Polisher



Order code:AM-040

ATR-E1 high frequency Electricity Knife



Order code:AM-044

LED vacuum cleaner



Order code:AM-049

Aluminum Oxi



0μm 90μm 30μm 460g

Order code:AM-041

Two pen air polisher



Order code:AM-045

Air polisher



Order code:AM-050

Teeth whitening machine

40W Teeth Whitening machine



Input voltage AC110V~240V,50/60Hz
The output voltage 18V
power 40W
Blue wavelength 460-490nm
Blue light density 200-300mW/cm²
Number of LED tubes 8
Package dimensions
Package weight 69x69x21cm
Packing quantity 1set/box
Whitening time Irradiate 2 to 3 times per
treatment course 10-15min each time

65W Teeth whitening accelerator



Product name:Teeth Whitening Machine
Light mode:blue;red; blue+red+purple;blue+purple
Board spectrum:450~510nm (blue light)
620-640nm (red light)
380-420nm (purple light)
Light source:4~5W / pc*8pcs (blue light)
2-3W / pc*4pcs (red light)
1W / pc*5pcs (purple light)
Voltage (to power supply):100-230VAC
Frequency:50-60Hz
Light Size:75*25mm
Operation Time:0-30mins
Size (L*W*H):75*75*150cm
Single package size: 96X44X19 cm
Single gross weight:9.0 kg

Wave Length	Light Source	LED Diode QTY(pcs)	Power	Function
460nm±20	Blue Light	8	40W	General Bleaching
625nm+20	Red Light	4	20W	Bleaching Relax
395nm±20	Purple Light	5	5W	H2O2/TiO2 Gel Bleaching
Total		17	65W	

40W Teeth Whitening machine




Input voltage AC110V~240V,50/60Hz
The output voltage 18V
power 40W
Blue wavelength 460-490nm
Blue light density 200-300mW/cm²
Number of LED tubes 8
Package dimensions
Package weight 69x69x21cm
Packing quantity 1set/box
Whitening time Irradiate 2 to 3 times per
treatment course 10-15min each time
LED Bulb:blue led 4pcs Red led 4pcs

60W Teeth whitening accelerator



Input voltage Ac 100V~240V
Work rate 60W
The output voltage DC18V
Blue light wavelength 460-490nm
Blue light density 300~400mw/cm²
Constant temperature 35~40°C
Number of LED tubes 12
Whitening time 10-15minutes for one time,
two or three times for a treatment.
Weight 27KG
Package dimensions 80x58*30cm
Package weight 30KG



37°C

Controllable temperature



Teeth whitening machine

46W three color whitening accelerator

54W three color whitening accelerator

18W three color whitening accelerator

SMD whitening accelerator



Desk version Chair version

Order code:AM-051



Desk version Chair version

Order code:AM-052



Desk version Chair version

Order code:AM-053



Light color blue
Power outputs up to 2500 mW/cm
Voltage 100/240 Volt, 50/60 Hz
Broad spectrum 100%, 80%, 60%, 40%
Light source 430nm ~ 490nm (Blue light)
Light timing 6 high power 3 watt LEDs *4pcs
Size 1~30 mins

Order code:AM-054

Three color whitening accelerator

36W whitening accelerator



Order code:AM-055

Light color blue, red, blue+red+purple,blue+purple
Voltage AC100-240V,50/60Hz
Led light 15pcs
Wavelength 430nm~490nm for Blue lights
620nm~640nm for RED lights
380nm~420nm for purple lights
Light source 4~5W/pc(Purple Light)*5pcs;
2~3W/pc(Red Light)* 2 pcs;
1W/pc(Purple Light)*8 pcs.
Light timing 1~30 min
Packing Size 69X48X21 cm
Gross Weight 9.5kg



Order code:AM-056

Light color Blue light/ Blue+Red+Purple light
Power supply AC100-240V,50/60Hz
Power 36 watt
Wavelength 420-480nm
Light size 18*74mm
Light intensity > 6000mw/cm
Light timing 0~20 min/0-40min
Power instabil <1%
Working met time output & continuous working
Coupling

Teeth whitening machine

3ml/5ml/10ml whitening gel



Order code:AM-057

desensitizing gel



Order code:AM-058

2.5ml/5ml Dual gel



Order code:AM-059

Gum protector



Order code:AM-060

Whitening kit



Order code:AM-061

Cheek retractor



Order code:AM-062

Whitening kit




Order code:AM-063

Whitening kit



Order code:AM-064

Whitening kit



Order code:AM-065

Oral brush up



Order code:AM-066

Dental bibs



Order code:AM-067

Cotton rolls



Order code:AM-068

3D Shade guide



Order code:AM-069

Disposable Saliva towel



Order code:AM-070

protect glasses



Order code:AM-071



Order code:AM-072

Vitamin E




Order code:AM-073

Root Apex Locator



Parameter:

Battery : 3.7V/750mAh
Power consumption ≤ 0.5W
Screen : 4.5" LCD
Weight : 350g
Dimensions : 82mm*88mm*112mm



Order code:AM-074

Root Apex Locator



Parameter:

Battery 3.7V/750mAh
Voltage 100V-240V 50Hz/60Hz
Dimensions 131 × 75 × 105mm (H)
Screen 3.2" LCD
Power source Electricity
Warranty One year



Order code:AM-075

Root Apex Locator



Parameter:

Battery:DC3.7V/750mAh
Current:170mA
Working system:Short-time operation
weight:350g
Dimensions : 84*88*112mm
Display : 4.5" LCD
Gross weight 0.7kgs
Package dimensi 17*14*9.5cm



Order code:AM-076

Root Apex Locator



Parameter:

Battery:ICR14500,3.7V 800mAh
Power Adapter:AC100-240V,50/60Hz
Input power:5VA
Screen size:4.0 inch
weight:350g




Order code:AM-077

Root Apex Locator



Parameter

Battery: 4.5V, 3 pcs of AAA batteries
Display screen: 3.7 inch
Main body size: 80*60*55mm



Order code: AM-078

AI Apex locator



Parameter

Battery: 3.7V/2000mAh
Adapter: Input: ~100V-240V 50Hz/60Hz 0.4A
Max Output: DC5V/1A 1.93
Consumption power: ≤0.5W 1.94
Screen: 3.8" LCD




Order code: AM-080

Refine A7 Apex locator



Parameter

Input (charge the battery): DC 5V 1A
AC Adapter: 100-240V AC 50Hz/60Hz
Battery: Li-ion Battery DC 3.6V/850 mAh
Working Voltage DC 3.6V
Working Current DC: 0.13-0.14A
Rated Power ≤0.5W
Measurement Voltage AC 200mV
Current AC 100μA
Screen 5.1" LCD (110 mm x 70 mm)



Order code: AM-079

RPEX6 Apex locator



Parameter

Battery : 3.7V/2000mAh
Adapter : ~100V-240V 50Hz/60Hz
Consumption power ≤0.5W
Screen : 4.5" LCD
Operation condition
Environment temperature : 0~40°C
Relative humidity : 10~85%RH.
Atmosphere pressure 60kPa~106kPa
Dimensions : 131mm x 75mm x 105mm



Order code: AM-081

Apex locator with endo motor



Parameter

Chuck type	Push button
Gear ratio	16:01
Max speed	600 rpm
Spray	no water
Max torque	4.0 N.M
Noise	≤55db



Order code:AM-082

C-SMART-I Endo motors



Parameter

Power Adaptor :Input:	100-240Vac, 50/60Hz 0.8A
Output:	DC10V/1.5A
Battery: Voltage:	DC7.4V
Capacity:	1000mAh
Rotation speed:	125-625rpm
Adjustable torque:	0.3-5.2N.cm
Volume(cm):	24x21x9.5
Weight(kg)/PCS:	1.8

Order code:AM-084

C-SMART-I PILOT Endo motor with Apex locator



Parameter

Adaptor Input:	AC100~240V, 50/60Hz
Output:	DC10V, 1.5A
Control unit's battery:	Li-ion battery(DC7.4v 2600mAh)
Handpiece's battery:	Li-ion battery(DC3.7v 1200mAh)
Rotation Speed range:	150~1000rpm
Torque:	0.6~3.9N.cm
Input power:	35VA
Volume(cm):	24x24x7.5
Weight(kg)/PCS:	1.5



Order code:AM-083

C-SMART-I PRO Endo motors



Parameter

Adapter Input:	AC100-240V, 50/60Hz
lithium battery:	DC7.4V, 2600mAh
Speed range:	150-650rpm
Torque range:	0.6-5.2N.cm
Protection against electrical shock:	Type B
Volume(cm):	23.5x20x14.5
Weight(kg)/PCS:	2.8



Order code:AM-085

Wireless Endomotor with LED Light



Parameter
Voltage:~AC100-240V , 50/60Hz
speed:16:1
Torque range:0.6-3.9N.cm (6-39mN.m)
Speed range:125rpm-625rpm
battery capacity:1200mAh
Illumination:5000lux
Operating mode:6 modes
Package dimensions:13.5*18*7.5
Package gross weight:0.5kgs



Order code:AM-086

Wireless Endomotor with LED Light



Parameter
Voltage:~AC100-240V , 50/60Hz
speed:16:1
Torque range:0.6-3.9N.cm (6-39mN.m)
Speed range:125rpm-625rpm
battery capacity:1200mAh
Illumination:5000lux
Operating mode:6 modes
Package dimensions:13.5*18*7.5
Package gross weight:0.5kgs



Order code:AM-087

2 in 1 Endomotor with apex locator



Features

1. With Multi-frequency length measuring technology and the function of root-canal length measuring and enlargement
2. Inside Apex locator Can be used alone for root measurement. Requency multiplexing technique and high precision measurement
3. 10 memory modes of adjustable speed,torque and angle. Clockwise rotation, negative rotation,reciprocating mode, root measurement,torque control and so on.
4. Continuous rotation mode Ajustable speed: 160-1000rpm
5. Auto reverse Angle adjustable
6. The endomotor is combined with apex locator.

Parameter

Lithium battery	3.7V/8000mAh
Torque range	0.5-3.0Ncm
Range of rotation	120-1000 rpm

Order code:AM-088

2 in 1 Endomotor with apex locator



Features

1. With Multi-frequency length measuring technology and the function of root-canal length measuring and enlargement
2. Inside Apex locator Can be used alone for root measurement. Requency multiplexing technique and high precision measurement
3. 10 memory modes of adjustable speed,torque and angle. Clockwise rotation, negative rotation,reciprocating mode, root measurement,torque control and so on.
4. Continuous rotation mode Ajustable speed: 160-1000rpm
5. Auto reverse Angle adjustable
6. The endomotor is combined with apex locator.

Parameter

Lithium battery	3.7V/8000mAh
Torque range	0.5-3.0Ncm
Range of rotation	120-1000 rpm

Order code:AM-089

16:1 Wireless Endomotor with LED Light



Parameter

Power Adapter:AC100-240V
Torque:0.6N.cm-4.0N.cm
Rotate speed:120rpm-700rpm
Round trip angle
30° -150°,150°-30°,90°-90°




Order code:AM-090

1:1 LED Wireless Endomotor with LED Light



Parameter

Power Adapter:AC100-240V
Torque:0.6N.cm-4.0N.cm
Rotate speed:120rpm-1000rpm
REC adjustable angle
30°,60°,90°,150°,180°, 210°,250°,370°
Rate:1:1



Order code:AM-091

2 in 1 16:1 Endomotor with apex locator



Parameter

Power Adapter:AC100-240V
Torque:0.6N.cm-4.0N.cm
Rotate speed:120rpm-650rpm
REC adjustable angle
30°,60°,90°,150°,180°, 210°,250°,370°




Order code:AM-092

2 in 1 1:1 Endomotor with apex locator



Parameter

Power Adapter:AC100-240V
Torque:0.6N.cm-4.0N.cm
Rotate speed:120rpm-1000rpm
REC adjustable angle
30°,60°,90°,150°,180°, 210°,250°,370°
Rate:1:1



Order code:AM-093

Prophy Pro



Technical:

Power Adapter;AC100-240V
Torque:>1.5N.cm
Speed range:500rpm-2500rpm
Lithium battery:3.7V/1900mAh
Screen:0.96 inch TFT color LCD

Ultra X cordless ultrasonic activator



Order code:AM-097

Dental Gutta Percha Obturation Pen




Order code:AM-095

Order code:AM-096

Technical :

- 1.Battery:DC3.7v,2800mAh
- 2.Heating time:2seconds
- 3.Working temperature:150 c/200 c
- 4.Adpter:Input 100v-240v AC ,Output DC5v,2.4v



Endo ultra activator



Features

- Frequency 40KHz
- One key work, indicator light prompt
- Head can rotate 360 degree
- Main unit only 120g weight;
- Strong working power,big capacity battery
- Safe enough for no tip break
- USB charging



Order code:AM-099

Gutta cutter



Order code:AM-100

Gutta cutter



Order code:AM-101

Pulp Tester



Order code:AM-102

Feature:

- 1.Preset speed mode(high-mid-low speed)
- 2.Agentle,pulsed stlmulus begins to Increase at a rate specifically designed for patient comfort.
- 3.If the patient indicates perception,simly release the button, the stimulus stops immediately,but the numbers remain frozen on the screen for about 3 minutes.
- 4.Specififfically designed for patient comfort.
- 5.Convenient to operate.
- 6.Maximum stimulus- - 80.between 0-40,the patient feels ache and anesthesla,alive teeth nerve.
- 7.Between 40 80, with above-mentioned reaction,part of teeth nerve dead;80,no above-mentioned reaction,dead teeth nerve.
- 8.The device turns itself offfff in three minutes after operation

Caries Detector



Order code:AM-103

Technical:

Adaptor working voltage:100-240Vac, 50/60Hz;
Lightwave Data: LED intensities:40-280mw
LED wavelength:380-420nm
Volume(cm):16X14X6.5
Weight/PCS:0.5 kg
Packing material:carton
QTY per master carton:6 pc/ctn



Fog mirror



Order code:AM-104



Base Light



Order code:AM-105

Feature:

- 1.Effiffifficient shade-taking without leaving the dental chair,good for dentist
- 2.Avoiding dentist choosing the wrong teeth color due to the unstable light
- 3.Three light setting:day light/room light/ ambient light,easily switch
- 4.12 pcs super LED light



5500K

3200K

3900K

C-FILL mini Endodontic obturation system




Technical:
Adapter: Input AC100~240V 50/60Hz
Output: DC 5V,1.5A
Battery: Rechargeable Li-ion battery
(DC3.7V 2000mAH)

Order code:AM-106

Obturation system



Technical:
Adapter: Input AC100~240V 50/60Hz
Output: DC 5V,1.5A
Battery: Rechargeable Li-ion battery
(DC3.7V 2000mAH)



Order code:AM-108

FI-G with FI-P Obturation system



Technical:
Adapter: Input AC100~240V 50/60Hz
Output: DC 5V,1.5A
Battery: Rechargeable Li-ion battery
(DC3.7V 2000mAH)

Order code:AM-107

Fast-Fill& Fast-Pro Obturation system



Technical:
Adapter: Input AC100~240V 50/60Hz
Output: DC 5V,1.5A
Battery: Rechargeable Li-ion battery
(DC3.7V 2000mAH)



Order code:AM-109

5S curing light



Order code:AM-110

Parameter:

Adaptor Input : ~110V-240V 50-60Hz
Adaptor output : DC4.2V 1A
Battery voltage : DC3.7V
Battery Volume : 1500mAh
wavelength : 420nm-480nm
Intensity : 1200-1500mw/cm2
Curing time : 5s-40s
Pacakge dimension : 24*17.5*10.5cm
Weight : 485g

5W 5s curing light



Order code:AM-111

Parameter:

Working mode : Strong、Pulse、渐亮
Curing time : 5s-40s
Curing depth : 5s/3mm
Intensity : $\geq 1600\text{mw/cm}^2$
Battery Volume : 2200mAh/3.7v
Input voltage : AC100-240V 50/60Hz
Package : 23*17.5*55cm
Gross weight : 0.57kg



Order code:AM-112

Parameter:

Adaptor input : ~110V-240V 50-60Hz
Adaptor output : DC4.2V 1A
Battery voltage : DC3.7V
Battery Volume : 2200mAh
wavelength : 420nm-480nm
Intensity : 1200-1600mw/cm2
Package dimension : 21.5*15*6cm
Weight : 0.7KG



Order code:AM-113

Parameter:

Adaptor input : ~110V-240V 50-60Hz
Adaptor output : DC4.2V 1A
Battery voltage : DC3.7V
Battery Volume : 2200mAh
wavelength : 420nm-480nm
Intensity : 1200-1600mw/cm2
Package dimension : 21.5*15*6cm
Weight : 0.65KG



Order code:AM-114

Parameter:

Adaptor input : ~110V-240V 50-60Hz
Adaptor output : DC4.2V 1A
Battery voltage : DC3.7V
Battery Volume : 2200mAh
Wavelength : 420nm-480nm
Intensity : 1200-1600mw/cm2
Package dimension : 21.5*15*6cm
Weight : 0.65KG



Order code:AM-115

Parameter:

Adaptor input : ~110V-240V 50-60Hz
Adaptor output : DC4.2V 1A
Battery voltage : DC3.7V
Battery Volume : 1500mAh
wavelength : 420nm-480nm
Intensity : 1200-2000mw/cm2
Package dimension : 21.5*15*6cm
Weight : 0.35KG



Order code:AM-116

Parameter:

Adaptor input : ~110V-240V 50-60Hz
Adaptor output : DC4.2V 1A
Battery voltage : DC3.7V
Battery Volume : 2200mAh
wavelength : 420nm-480nm
Curing time:5s-40s
Package : 20.5*13.5*6.8cm
Weight : 0.341kg



Order code:AM-117

Parameter:

Watts : 7W
Adaptor input : ~110V-240V 50-60Hz
Adaptor output : DC4.2V 1A
Battery voltage : DC3.7V
Battery Volume : 2200mAh
Wavelength : 440nm-490nm
Intensity : 1200-1500mw/cm2
Curing time : 5S-40s
Package : 24*7.5*6cm
Weight : 0.545kgs

Maxcure 9 1s curing light



Order code:AM-118


Model: MaxCure 9
Voltage: 100-240VAC 50Hz/60Hz 0.4A
Output: DC 5V/1A
Wavelength: 380nm ~ 515nm
Light intensity: 1000mW/cm² ~ 2500mW/cm²
Battery: lithium battery 2600mAh
Package size : 20X15X10 cm
Gross weight : 2.0 kg

VRN VAFU 1S curing light



Order code:AM-119

Model: VRN VAFU
Voltage: 100-240VAC 50Hz/60Hz 0.4A
Output: DC 5V/1A
Wavelength: 380nm ~ 515nm
Light intensity: 1600mW/cm² ~ 3200mW/cm²
Battery: lithium battery 2600mAh
Package size : 20X15X10 cm
Gross weight : 2.0 kg



10W 1s Curing light




Order code:AM-120

Parameter
Watts:10W
Voltage: 3.7-4.2V
Mode: strong light, pulse, fade
Curing time: 1S
Curing depth: 1s/3mm
Light intensity: >2000mw/cm²
Battery capacity: 2200mAh
Product weight: 0.75kg
Packing size: 26*16.5*6.5cm



10W 1s Curing light



Order code:AM-121

Working mode	Light intensity mW/cm ²	
	2500	1200
Normal	1s,2s,3s	3s,5s,10s , 15s
Pulse		3s,5s,10s
Ramp		5s,10s,15s,20s
Power dapter	AC100-240V	
Li Battery	18650/2600mAh	
4 LED(Blue)	10W	
Wave length	420-515nm	



10W 1s Curing light



Order code:AM-122

Dental chair curing light



Order code:AM-124

Working mode	Light intensity mW/cm ²	
	2500	1200
Normal	1s,2s,3s	3s,5s,10s , 15s
Pulse		3s,5s,10s
Ramp		5s,10s,15s,20s
Power dapter	AC100-240V	
Li Battery	18650/2600mAh	
4 LED(Blue)	10W	
Wave length	420-515nm	



Light guide



Order code:AM-125

Wide spectrum 1s Curing light



Order code:AM-123


Curing light with Caries Detector



Order code:AM-125


Parameter

Battery 2200MA/H3.
Adapter AC100-240V
Battery DC5V 1A
Wave length 385-515nm
Atmosphere 70kPa-
Package size 268.8*26.1*3
Net weight 259g
Light intensity : 1600-2700mw/cm2



Parameter

Output power : 5W
Charging voltage : AC100~240V,50/60Hz
Wave length : 420-480nm
Light intensity : 800-1800mw/cm2
Lithium battery : 3.7V 2200mAh
Single Dimension : 24.5*177*13.5cm
Weight:0.6kgs



Blue curing Purple detector

Dental Oral Microscope



Order code:AM-126

Parameter

Eyepieces : Binoculars 10.5X 15B(Aperture: 15mm)
 PD adjustable range : 55mm-80mm
 Diopter Range : $\pm 7D$
 Big objective lens : F=250mm
 Object diameter : LED light source>45,000LX
 View diameter : 10mm-100mm
 Spot diameter : 90mm
 Balance arm Down: 240mm,Middle: 450mm,Up: 230mm
 Maximum resolution : 1920*1080P
 Total magnification ratio : 3.35-22.5TIMES , 1600W
 Camera pixel : 1600W
 Working distance : about 250mm
 Microscope light source : LED Balance the light source
 Camera shooting : Extraposition
 Light color: white and orange
 Connection:computer, dental camera, monitor, TV, etc



OMS2000 Dental Oral Microscope



Order code:AM-127



LED Dental Implanting Lamp



Ceiling type double head



Order code:AM-128

Input voltage:DC12-24V
Power:30W(One head)
Illumination:30000-90000lx
Color Temperature:4000-5500k(Adjustable)
Spot Size:Φ120mm
LED Bulb:26pcs

Ceiling type single head



Order code:AM-129

Floor standing type



Order code:AM-130

Built-in HD camera

Pixel:2,000,000
Magnification:20
Optical Zoom:10
Image resolution:1080P
Storage:10GB
Wireless net:2.4G/5.8G
Voice control:Support
Remote video:Support
Online update:Support
LAN projection:Support



Ceiling type single head



Order code:AM-131

Input voltage:DC12-24V
Power:30W(One head)
Illumination:40000-90000lx
Color Temperature:4000-5500k(Adjustable)
Spot Size:Φ120mm
LED Bulb:32pcs

Ceiling type double head



Order code:AM-132



Dental unit type



Order code:AM-133

LED Dental Implanting Lamp



Order code:AM-136

Input voltage:DC12-24V

Power:12W

Illumination:8000-30000lx

Color Temperature:4000-5500k(Adjustable)

Spot Size:80*160mm(Distance 700mm)

LED Bulb:8pcs

Dental unit type



Order code:AM-134

Input voltage:DC12-24V

Power:25W

Illumination:8000-80000lx

Color Temperature:4000-5500k(Adjustable)

Spot Size:Φ120mm

LED Bulb:18pcs



Order code:AM-137

Input voltage:DC12-24V

Power:9W

Illumination:8000-30000lx

Color Temperature:4000-5500k(Adjustable)

Spot Size:80*160mm(Distance 700mm)

LED Bulb:6pcs

Ceiling type



Order code:AM-135



Order code:AM-138

Input voltage:DC12-24V

Power:5W

Illumination:8000-30000lx

Color Temperature:4000-5500k(Adjustable)

Spot Size:80*160mm(Distance 700mm)

LED Bulb:2pcs

LED Dental Implanting Lamp



Order code:AM-139



Parameter

Input voltage:DC12-24V
Power:9W
Illumination:8000-80000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:80*160mm
(Distance 700mm)
LED Bulb:6pcs



Order code:AM-140



Parameter

Input voltage:DC12-24V
Power:15W
Illumination:8000-80000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:80*160mm
(Distance 700mm)
LED Bulb:8pcs



Order code:AM-141



Order code:AM-142



Order code:AM-143



Parameter

Input voltage:DC12-24V
Power:9W
Illumination:8000-80000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:80*160mm
(Distance 700mm)
LED Bulb:6pcs

Order code:AM-144



LED Dental Implanting Lamp



Order code:AM-145

Parameter

Input voltage:DC12-24V
Power:5W
Illumination:8000-30000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:80*160mm
(distance700mm)
LED Bulb:2pcs



Order code:AM-146

Parameter

Input voltage:DC12-24V
Power:25W
Illumination:30000-80000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:Φ120mm
LED Bulb:20pcs



Order code:AM-147

Parameter

Input voltage:DC12-24V
Power:15W
Illumination:8000-30000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:80*160mm
(distance700mm)
LED Bulb:8pcs



Parameter

Input voltage:DC12-24V
Power:15W
Illumination:8000-30000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:80*160mm
(distance700mm)
LED Bulb:6pcs

LED Dental Implanting Lamp

Double color implanting lamp



Order code:AM-149

Parameter

Input voltage:DC12-24V
Power:28W
Ilumination:8000-23000lx
Color Temperature:4000-5500C
(Adjustable)
Spot Size::80*160mm
(distance700mm)
LED Bulb:9*2pcs

Double color implanting lamp



Order code:AM-150

Parameter

Input voltage:DC12-24V
Power:18W
Ilumination:8000-20000lx
Color Temperature:4000-5500l
(Adjustable)
Spot Size::80*160mm
(distance700mm)
LED Bulb:6x2pcs

Double color implanting lamp



Order code:AM-151

Parameter

Input voltage:DC12-24V
Power:12W
Ilumination:8000-30000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:Φ120mm
(distance700mm)
LED Bulb:4x2pcs

Dental Implanting Lamp



Order code:AM-152

Parameter

Input voltage:DC12-24V
Power:9W
Ilumination:8000-30000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:Φ120mm
(distance700mm)
LED Bulb:6pcs

Dental Implanting Lamp



Order code:AM-153

Parameter

Input voltage:DC12-24V
Power6W
Ilumination:8000-25000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:Φ120mm
(distance700mm)
LED Bulb:4pcs

Double color implanting lamp



Order code:AM-154

Parameter

Input voltage:DC12-24V
Power:12W
Ilumination:8000-30000lx
Color Temperature:4000-5500k
(Adjustable)
Spot Size:Φ120mm
(distance700mm)
LED Bulb:1pcs

2.5X 3.5X dental loupes



Order code:AM-155



Order code:AM-156

2.5X 3.5X dental loupes



Order code:AM-157



Order code:AM-158

4X dental loupes



Order code:AM-159



Order code:AM-160

5X dental loupes



Order code:AM-161




Order code:AM-162

6X dental loupes



67 Order code:AM-163



Order code:AM-164

8X dental loupes



Order code:AM-165



Order code:AM-166

Dental chair



Order code:AM-167

Dental chair



Order code:AM-169

Dental chair



Order code:AM-168

Dental chair



Order code:AM-170

Dental chair



Order code:AM-171

Dental chair




Order code:AM-172

Dental chair



Order code:AM-173

Children Dental chair



Order code:AM-174

Portable dental unit



Order code:AM-175

Standard Accessories

Power:550W
Built- In air compressor...1套/set
Saliva ejector...1套/set
3- way air/ water syringe...1套/set
High speed air turbine Handpiece tube...1套/set
Low speed air turbine Handpiece tube...1套/set
Foot pedal...1套/set
600ML Clean water bottle...1套/set
600ML 3- way syringe bottle...1套/set
600ML Waste water bottle...1套/set
N.W./ G.W.:24Kg/28Kg
Dimension:L43*W34*H62cm
Package size:L51*W40*H69cm



Order code:AM-176

Standard Accessories

Power:550W
Built-In air compressor...1套/set
Saliva ejector...1套/set
3- way air/ water syringe...1套/set
High speed air turbine Handpiece tube...1套/set
Low speed air turbine Handpiece tube...1套/set
Foot pedal...1套/set
600ML Clean water bottle 600ML...1套
N.W./ G.W.:22Kg/24Kg
Dimension:L50*W27*H47cm
Package size :
L53*W34*H52cm



Order code:AM-177

Standard Accessories

Power:550W
Built-In air compressor...1套/set
Saliva ejector...1套/set
3- way air/ water syringe...1套/set
High speed air turbine Handpiece tube...1套/set
Low speed air turbine Handpiece tube...1套/set
Foot pedal...1套/set
600ML Clean water bottle 600ML...1套/set
N.W./ G.W.:22Kg/24Kg
Dimension:L50*W27*H47cm
Package size :
L53*W34*H52cm



Order code:AM-178

Standard Accessories

Power:600W
Pressure:110v-220V/50Hz-60Hz
Speed R.P.M:1400r/min
Air Flow:65L/min
Pressure:0.6Mpa
Noise:39dB(A)
Standard Accessories:
High speed air turbine Handpiece tube...1set
Low speed air turbine Handpiece tube...1set
3- way air/ water syringe...1set
Water bottle...1set
Drain bottle...1pcs
Suction System...1set
Foot Control...1pcs



Order code:AM-179

Standard Accessories

Power:600W
Pressure:110v-220V/50Hz-60Hz
Speed R.P.M:1400r/min
Air Flow:65L/min
Pressure:0.6Mpa
Noise:39dB(A)
Standard Accessories:
High speed air turbine Handpiece tube...1set
Low speed air turbine Handpiece tube...1set
3- way air/ water syringe...1set
Water bottle...1set
Drain bottle...1pcs
Suction System...1set
Foot Control...1pcs



Order code:AM-180

Standard Accessories

Power:600W
Pressure:110v-220V/50Hz-60Hz
Speed R.P.M:1400r/min
Air Flow:130L/min
Pressure:0.8Mpa
Noise:40dB(A)
Standard Accessories:
High speed air turbine Handpiece tube...1set
Low speed air turbine Handpiece tube...1set
3- way air/ water syringe...1set
Water bottle...1set Drain bottle...1pcs
Suction System...1set Foot Control...1pcs
Built-in Ultrasonic scaler...1pcs
Built-in LED Curing light...1pcs
6L Air tank...1pcs



Order code:AM-181



Order code:AM-182



Order code:AM-183



Order code:AM-184



Order code:AM-185



Order code:AM-186



Order code:AM-187



Order code:AM-188

Portable dental chair



Order code.:AM-233



Order code.:AM-234



Order code.:AM-235



Order code.:AM-236



Order code.:AM-237

CLASSIC FOUR COLORS ARE RECOMMENDED FOR YOU



Air compressor

One to one Air compressor



Features

Voltage: 220VAC 50Hz
Electric current: 2.4A
Power: 545W
Frequency: 70L/min
Start pressure: 0.5Mpa
Capacity pressure: 0.8Mpa
Air tank volume: 32L
Noise: 55-62dB
G.W.: 27kg
Packing size: 46x46x68cm (carton)

Order code.:AM-189

One to Two Air compressor



Features

Voltage: 220VAC 50Hz
Electric current: 3.7A
Power: 840W
Frequency: 105L/min
Start pressure: 0.5Mpa
Capacity pressure: 0.8Mpa
Air tank volume: 32L
Noise: 58-64dB
G.W.: 29Kg
Packing:Carton box

Order code.:AM-190

One to Two Air compressor



Features

Voltage: 220VAC 50Hz
Electric current: 3.7A
Power: 840W
Frequency: 105L/min
Start pressure: 0.5Mpa
Capacity pressure: 0.8Mpa
Air tank volume:40L
Noise: 58-64dB
G.W.: 32Kg/33Kg
Packing:Carton box
Packing size: 46x46x74cm(38L)

Order code.:AM-191



Features

Voltage: 220VAC 50Hz
Electric current: 5A
Power: 1 100W
Frequency: 1 40L/min
Pressure: 0.5Mpa
Capacity pressure: 0.8Mpa
Air tank volume: 45L
Noise:63-65dB
G.W: 38kg
Packing size: 53x53x77cm
(carton packing)

Order code.:AM-192

One to Four Air compressor



Features

Voltage: 220VAC 50Hz
Electric current: 7.4A
Power: 1680W
Frequency: 210L/min
Start pressure: 0.5Mpa
Capacity pressure: 0.8Mpa
Air tank volume: 65L
Noise: 70dB
G.W.: 61kg
Packing size: 60x60x82cm (carton)

Order code.:AM-193



Features

Electric Current:4.8A
Frequency:1090W
Power:140L/min
Start Pressure:0.5Mpa
Capacity Discharge:0.8Mpa
Air Container Volume:45L
Noise:62-66dB
G.W.:52KG
Packing Size:72cm x 44cm x 77cm (Wooden Box)

Order code.:AM-194



Features

Electric current: 7.4A
Power: 1680W
Frequency: 210L/min
Start pressure: 0.5Mpa
Capacity pressure: 0.8Mpa
Air tank volume: 60L
Noise: 65-68dB
G.W: 57kg
Packing size: 86x45x74cm(Carton box)

Order code.:AM-195

750L Dental suction machine



Features

Voltage(V): 220V 50-60Hz
Power(W):800W/850W
Rated Current(A):4.8/5.2
Max Vacuum(Kpa):≤-15
Air Flow (L/min)≥500
Product weight (kg): 26
Package Weight (kg)) : 37
Package size (cm):53*37*65
Noise(dB):67

77 Order code.:AM-196

1100L Dental suction machine



Features

Voltage(V): 220V 50-60Hz
Power(W):1200W/1400W
Rated Current(A):7.3/8.3
Max Vacuum(Kpa):≤-15
Air Flow (L/min)≥800
Product weight (kg): 34.5
Package Weight (kg)) : 46
Package size (cm):55*39*72
Noise(dB):69

Order code.:AM-197

1300L Dental suction machine



Features

Voltage(V): 220V 50-60Hz
Power(W):1300W/1500W
Rated Current(A):3.5/3.8
Max Vacuum(Kpa):≤-30
Air Flow (L/min)≥800
Product weight (kg): 41
Package Weight (kg)) :59.5
Package size (cm):56*48*79
Noise(dB):69

Order code.:AM-198

Dental suction machine



Features

Voltage(V): 220V 50-60Hz
Power(W):600W/620W
Rated Current(A):3.7/4.5
Max Vacuum(Kpa):≤-12
Air Flow (L/min)≥300
Product weight (kg): 32
Package Weight (kg)) :35
Package size (cm):45*43*96
Noise(dB):64

Order code.:AM-199

Dental cabinet



Table Color choice



Six drawers cabinet
Order code.:AM-200



Lock bar Selective Color



Stainless Six drawers cabinet
Order code.:AM-201



Table Color choice



Six drawers cabinet
Order code.:AM-202



Table Color choice



Five drawers cabinet
Order code.:AM-203

Dental cabinet



Lock bar Selective Color



Lock bar Selective Color



Six drawers cabinet

Order code.:AM-204

Six drawers cabinet

Order code.:AM-205



Cabinet body color options



Six drawers cabinet

Order code.:AM-205C



Lock bar Selective Color



Order code.:AM-205P



Table color



Order code.:AM-205F



Order code.:AM-205S



Order code.:AM-206



Order code.:AM-207

Dental cabinet



Order code.:AM-208



Order code.:AM-208A



Order code.:AM-208B

Three-layer steel cart



Order code.:AM-209



Order code.:AM-210



Order code.:AM-211



Order code.:AM-2112

Dentist chair



Order code.:AM-212



Order code.:AM-213



Order code.:AM-214



Order code.:AM-215



Order code.:AM-216



Order code.:AM-217



Order code.:AM-218



Order code.:AM-219



Order code.:AM-220



Order code.:AM-221



Order code.:AM-222



Order code.:AM-223



Order code.:AM-224



Order code.:AM-225



Order code.:AM-226



Order code.:AM-227



Order code.:AM-228



Order code.:AM-229



Order code.:AM-230



Order code.:AM-231



Order code.:AM-232

Micromotor



Type : Brushless Micromotor
Input voltage : 110/240V
Rotation : 60.000rpm
Max torque : 7.5N.cm
Weight : 4kg
Dimension : 34.5*26*20cm

Brushless Micromotor S200
Order code::AM-238



Micromotor 204
Order code::AM-241



Micromotor(204+102L)
Order code::AM-242



Micromotor(90+108E)
Order code::AM-243



Micromotor(90+102)
Order code::AM-240



70,000 rpm Micromotor
Order code::AM-244

Power : 150W
Max Torque : 800GF.Cm
Rotated speed : 0-35,000 RPM
Voltage:AC110V/220V 50/60HZ
(110V need to add a adapter)
Bur dimension:φ2.35mm
Power Supply voltage:110V-220V
Cycle:50/60 Hz
Electric Power:50VA
R.P.M:Handpiece Max

Power Supply voltage:110V-220V
Cycle:50/60 Hz
Electric Power:50VA
R.P.M:Handpiece Max
Weight:2.9kg
Dimension:157(W) X 202(D) X 250(H)mm
Nonload:1,000~70,000rpm
Max Torque:8.5Ncm
Bur size:2.35mm
Electric Power:230W
Net Weight:238g
Dimension:159(L)mm X 27(D)mm
Gross weight:2.5kg



Micromotor
Order code::AM-245



Order code::AM-246



70,000 rpm Micromotor
Order code::AM-247



Order code::AM-248



60,000rpm Micromotor
Order code::AM-249



Micromotor
Order code::AM-250



Micromotor
Order code::AM-251



204 micromotor handpiece
Order code::AM-252



90 micromotor handpiece
Order code::AM-253



SDE-M33Es micromotor handpiece
Order code::AM-254



SDE-H37LN
Order code::AM-255



Handpiece cable
Order code::AM-256



Foot Pedal
Order code::AM-257

C-Pumas master Brushless Electric Micromotor



Input power:AC100-240V, 50/60Hz
Display:5 inch LCD screen, 800*480px
Maximum speed:40,000rpm
Motor outer diameter:20.6mm
LED brightness:> 25,000lux
Speed:10,000~200,000rpm

Order code.:AM-258

Electric Brushless Micromotor



Power supply:AC110-220V 50/60HZ
Power:120w
Speed range:
2000rpm40000rpm (1:1)
100rpm-2500rpm(16:1)
10000rpm-200000rpm(1:5)
LED maximum:> 25000lux
Torque output:3.5N.cm
Volume:29X23X11cm

Order code.:AM-259

Elite Implanting machine



Power supply:AC110-220V 50/60HZ
Watts:210W
Motor Rotation Speed Range:1000rpm 30000rpm
Torque Range:5.0-55N.cm (20:1)
Stored Programs:10 types

Order code.:AM-260

C-SAILOR Dental implantation system



Power supply:AC230V 50/60Hz
Maximum output power:120W
Adjustable torque ,range:
5.0-55N.cm(20:1)
BLDC Speed:300rpm-50000rpm
Infusion flow(maximum):150ml/min
Motor:surgical brushless motor
Volume(cm):40X33X26.5
Weight(kg):8.2kg
Packing material: Aluminum box

Order code.:AM-261

C sailor Pro Dental implantation system



Power supply : AC220-230V 50/60Hz
Motor speed : 300--40000rpm
Adjustable torque :
5Ncm--80Ncm
Water supply volume :
4 blocks adjustable. 150ml/min
Light source illumination :
30000lux
Motor : surgical brushless motor
Volume(cm) : 42x38x26cm
Weight(kg) : 5.5kg

Order code.:AM-262

Implanting system



Voltage:AC 220V-240V 50-60HZ
Voltage :± 10Fluctuation
Frequency:50 Hz
Max. :340VAPower Consumption
Rated Voltage Torque Range:0~80N.cm
Rated Voltage Speed Range:
300~40,000 rpm
Cooling Water Flow:60/70/80 ml/min
Weight(Including Motor)(kg): 5.4
Dimension(mm):300×280×128

Order code.:AM-263

Implanting system



Voltage:AC 220V-240V 50-60HZ
Voltage :± 10Fluctuation
Frequency:50 Hz
Max. :340VAPower Consumption
Rated Voltage Torque Range:0~80N.cm
Rated Voltage Speed Range:
300~40,000 rpm
Cooling Water Flow:60/70/80 ml/min
Weight(Including Motor)(kg): 5.4
Dimension(mm):300×280×128

Order code.:AM-264

Ultrasonic Dental Piezo Bone Surgery Unit



Power-supply voltage:
100~240V 50/60Hz
Watts:150VA
Fuses:2×F10A 250V
Working frequency:
24kHz ~ 36kHz
Water Flow:60~100ml/min

Order code.:AM-265

MaxSurgery III UltrasonicBoneCutter



Mechanical Size: 290mm*270mm*125mm
Main Unit Weight: 3kg
Power Input: 100-240V A.C. 50Hz/60Hz
Main Unit Input: 120VA
Frequency: 26±3kHz
Water Flow Rate: 30~125ml/m
Operating Model: B, P,E, C
Control Model: Touch Control

Order code.:AM-266

ORAL CAMERA

Omega pixels CCD Intraoral Camera(USB+VGA)



Technical Parameter
 Display:17 inches
 Resolution:1920*1280
 Light:10 white LED light
 Pixel:10 million
 Handle cable:3M spiral line
 or 1.8M straight line
 Power supply:110-220V 12v/3A
 Interface:USB
 Power consumption:<30W
 Storage:U disk storage
 G.W:4.5KGS
 Packing:410*110*450mm

Order code.:AM-267

17" multi-media intraoral camera



Technical Parameter
 Display:17 inches
 Resolution:1920*1280
 Light:10 white LED light
 Pixel:10 million
 Handle cable:3M spiral line or
 1.8M straight line
 Power supply:110-220V 12v/3A
 Interface:USB
 Power consumption:<30W
 Storage:U disk storage
 G.W:4.5KGS
 Packing:410*110*450mm

Order code.:AM-268

17" multi-media intraoral camera



Technical Parameter
 Resolution:1280*1204
 Light:10 LCD cold light sources
 Pixel:10 million
 Handle cable:3M spiral line or
 1.8M straight line
 Power supply:110-220V 12v/3A
 Interface:USB
 Power consumption:<30W
 Storage:U disk storage
 G.W:4.5KGS
 Packing:410*110*450mm

Order code.:AM-269

17" multi-media intraoral camera



Technical Parameter
 Display:17 inches
 Resolution:1920*1280
 Light:10 white LED light
 Pixel:10 million
 Handle cable:3M spiral line or 1.8M straight line
 Power supply:110-220V 12v/3A
 Interface:USB
 Power consumption:<30W
 Storage:U disk storage
 G.W:4.5KGS
 Packing:410*110*450mm

Order code.:AM-270

Monitor screen



Monitor Bracket



WIFI oral camera



Pixel 500W
 Light freedom 6 white cold light sources
 Dimming 5-level brightness adjustment
 Focus range 10mm-50mm
 Connect 1 piece
 Available distance ≤30m
 Battery specifical 3.7V 1000mAh
 Charging time 2-3 hours
 Sustainable use ≥3h
 Connection methcWireless Wififi



Order code:AM-271



Order code:AM-272



Order code:AM-273



Order code:AM-274



Order code.:AM-275

ORAL CAMERA

Wireless Dental Endoscope Intraoral Camera



Order code:AM-276

Technical Parameter
Image sensor: 1/4 sony CCD
Image resolution: 2.0 mega pixels
signal system : NTSC/PAL
Lamp : 6pcs white LED (5600K)
Accumulation point : 3mm--50mm
Photographing visual angle: 105°
Lead : length 2.5m Connect port : 6 pin hole
Channel selection: 4 channel
Frequency: 5.8GHz
wireless distance : Over 10m
Power : DC5V±0.5 (1.5A)
Battery : 350mA lithium battery
(cold work about 35 mins continuously)

5.0mega pixels CCD Intraoral Camera(USB+VGA)



Order code:AM-277

Technical Parameter
Product name Split oral observation instrument
Image sensor: CMOS 1/4" HD sensor
System solution: embedded operating system
Display screen: single screen / four screens / sixteen screens
System operation mode: infrared remote control/key control
Backup method: automatic backup to U disk
Storage capacity: maximum support 32G storage
Picture storage format: JPG format
Hand shooting focus range: 5MM-50MM
WIF transmission: Direct wireless connection with computer
Package size : 12X22X22 cm

5.0mega pixels CCD Intraoral Camera(USB+VGA)



Order code:AM-278

Technical Parameter
Image Sensor:1/4" HD CMOS
Active Pixels:1280 x 1024 (1.3MP)
Resolution:1280 x 720 (720p)
FrameRate:30fps@720p
Focus Range:5mm~35mm
Angle of view:80°
Distortion: < 5%
Lighting:6 LEDs
Signal Output:USB2.0
Length of Cable:2m
Driver Type:UVC/WDM
TWIN:Yes
Operation System:
Windows2000/xp/win7
/win8/win8.1/win10 (32bit&64bit)

SD card video dental camera



Order code:AM-279

Technical Parameter
Image Sensor:1/4" CMOS sensor
Image Resolution:5.0 mega pixels
Output jack:Video
Lamp:6 pcs white LED(5600K)
Accumulation point:3~50mm
Photographing visual angle:105 degree
Signal format:NTS/PAL(opti)



Order code:AM-280

Dental intraoral camera

Technical Parameter
Camera diameter:8mm
Camera sensor:CMOS
Pixel resolution:1920 x 1080
Focal length:1-2cm
Waterproof rate:IP67
Battery capacity:1000mAh
Compatible system:
Android 4.2+/ iso 8.0+
Working time2 hours after
charging with full battery

Wifi Intraoral Camera(USB)



Order code:AM-281

Support for Android / Apple / Windows
Cellphone Pad Computer
Sensor: 1/2.7" CMOS
Resolution: 720P
Battery : 750 mA
Lithium battery
Light: Eight white LEDs

Touch screen X Ray



Order code:AM-282

Technical Parameter

Output power: Lithium battery 14.8VDC, 10A
Tube voltage/current: 60KV/1mAh (built-in adjustable)
Working frequency: 20KHz
Exposure time: 0.2-2.0S
Focus standard value: 0.3mm
Total capacity: 1.75mmAL
Display mode: LCD touch screen
Product size:19*11*17.5cm
Net weight: 2kg
Packing size:36*24.5*18.5cm
Gross weight:3.5kg

Handy X Ray



Order code:AM-284

Technical Parameter

Tube voltage/current: 70KV/2mAh (built-in adjustable)
Working frequency: 20KHz
Exposure time: 0.04-2.0S
Spot:0.4mm
Total capacity: 1.75mmAL
Product size:19*11*17.5cm
Net weight: 2kg
Packing size:36*27*20cm
Gross weight:4.1kg



X Ray



Order code:AM-283

Technical Parameter

Output power: Lithium battery 14.8VDC, 5000mA
Tube voltage/current: 60KV/1.5mAh (built-in adjustable)
Working frequency: 20KHz
Exposure time: 0.01-2.5S
Focus standard value: 0.3mm
Total capacity: 1mAL
Spot:0.4mm
Product size:19*11*17.5cm
Net weight: 2kg
Packing size:32*26*27cm
Gross weight:4.1kg

Refine X ray



Order code:AM-285

Technical Parameter

Tube Voltage : 60KV
Highest Voltage: 70KV
Tube Current: 2mA
Exposure Time: 0.01-2s
Focus: 0.4mm
Anode Target Angle: 12°
Total Filtration: 1.0mmAl
Battery Capacity: 950mAh
Product size:32*26*27cm
Net weight: 4.1kg
Packing size:32*26*27cm



Hyperlight Portable Touch screen X Ray



Order code:AM-286

Technical Parameter
 Charger power supply:~100-240V 50/60Hz
 Power supply:Lithium ion battery: 14.8V, 2500mAh
 Charger power output:24V -1.5A
 X-ray tube model:KL11-0.4-70
 Tube voltage:65kV($\pm 10\%$)
 Tube Current:2.5mA ($\pm 20\%$)
 Exposure time adjustment range:0.02-2S
 Nominal power:162.5W nominal at 65kV, 2.5mA
 Minimum source to skin distance:
 > 20cm (from focal spot to cone tip)



Order code:AM-287



Order code:AM-288



Order code:AM-289



Order code:AM-290



Order code:AM-291



Order code:AM-292



Order code:AM-293



Order code:AM-294



Order code:AM-295



Order code:AM-296



Order code:AM-297



Order code:AM-298



Order code:AM-299



Order code:AM-300




Order code:AM-301

Handy X ray sensor

Handy



Order code:AM-301



Specification	HDR-500	HDR-600
Chip Type	CMOS APS	
Fiber Optical Plate	Yes	
Scintillator	GOS/CSI	
Dimension	27.5x38.5	32.3x44.3
Active Area(mm)	22.5X30	27X36
Pixel Matrix	18.5	
Active Pixels	1.9M	2.7M
Operation System	Windows2000/xp/win7	
	win8/2in8.1/win10	
	(32bit&65bit)	

Ateco X ray sensor



Order code:AM-302

Item code:AT301

Active area:22.5 x 30mm

Sensing pixels:1.9Mega pixels

Launched software APP:iPhone & iPad /

Android system will be available



Refine X Ray sensor



Order code:AM-303

Specification	R1	R2
Detector Technology	APS CMOS	APS CMOS
Scintillator	CSL	CSL
Theoretical Resolution	25lp/mm	25lp/mm
True Resolution	20lp/mm	20lp/mm
Pixel Matrix	1000*1500	1300*1800
Active Area(mm)	20*30	26*36
Mechanical Size(mm)	25.4*36.8*4.4	30.4*41.9*4.4
Data interface	USB 2.0	USB 2.0
AD Conversion(bit)	16	16
Protection	IP68	IP68



Korea Vatech X ray sensor



Order code:AM-304



Detector:CMOS

Pixel Size:14.8μm

Theoretical Resolution:33.78lp/mm

Dynamic Range:12 bit

Dimensions(WxLXT):40.8 x 30.6 x 7.9 mm (1.5)

Active Area(WxL):33.00 x 23.98 mm (1.5)

Exterior Material:Silicone (Bio-compatibility)

Total Cable Length:2.7m


X Ray system

Nanopix X ray sensor



Order code:AM-305

Specification	Nanopix1	Nanopix2
Chip Type	CMOS APS	
Fiber Optical Plate	Yes	
Theoretical resolution	25lp/mm	
Scintillator	GOS/CSI	
Dimension	27.5x38.5	32.3x44.3
Active Area(mm)	20X30	26X36
Mechanical size(mm)	25.4*36.8*4.4	30.4*41.9*4.4
AD Conversion(bit)	16	
Operation System	Windows2000/xp/win7	
	win8/2in8.1/win10	
	(32bit&65bit)	



Handy X ray Scanner



Order code:AM-307



Technical Parameter
Pixel size:30*30μm
Scan Time:3s
Image Plate:Size 0/1/2/3
Driver:Twain
Weight:~1.4kg
Dimension (W*H*D):21.8*8*8cm
Operation system:
Windows7\8\10, Mac, Linux

Fussen X ray Scanner




Supported Photophor Size:Size 0: 22X31mm, Size 1: 24X 40mm, Size 2: 31 X 41mm, Size3: 27 X 54mm
Pixel Size:35um
USB Interface:USB2.0, USB3.0
Photophor Records:Automatic Clearing
Gray Scale:16bit
Theoretcal Resolution:130mm
Power Supply:100-240V(AC),50/60Hz,1.5A
Product Size:260×167×325mm
Weight:5.2kg
Packing Size:374×320×336mm
Software:Windows XP, Windows 7, Windows 8, Windows 10



Order code:AM-306

Panda P2 3D scanner



Scammomg Technology:Video Type
Overall Size(LXWXH):216X40X36mm
Tip Size(LXWXH):83.4X19.6X14.6mm
Data Format STL,PLY,PTY
Digital Interface:USB 3.0
Field of View:18x16mm
Depth of Field:20mm
Accuracy:10um
Weight:246g

Order code:AM-308

Steam Sterilizers

Digital Steam Sterilizers



16L/18L/23L

Order code:AM-309

LCD Steam Sterilizers



18L/23L

Order code:AM-310

LED Steam Sterilizers



17L/23L

Order code:AM-311

LCD Steam Sterilizers



17L/23L

Order code:AM-312

Simple Steam Sterilizers



18L

Order code:AM-313

Simple Steam Sterilizers



18L

Order code:AM-313

Simple Steam Sterilizers



18L

Order code:AM-314

COSOMO Steam Sterilizers



18L

Order code:AM-315

Disinfection cabinet



Order code:AM-316

Disinfection cabinet



Order code:AM-317

0.8L Ultrasonic Cleaner



Order code:AM-318

2.5L Ultrasonic Cleaner



Order code:AM-319

2.5L Ultrasonic Cleaner



Order code:AM-320

2.5L Ultrasonic Cleaner



Order code:AM-321

1.3L Ultrasonic Cleaner



Order code:AM-322

0.75L Ultrasonic Cleaner



Order code:AM-323

0.6L Ultrasonic Cleaner



Order code:AM-324

0.6L Ultrasonic Cleaner



Order code:AM-325

5L Ultrasonic Cleaner




Order code:AM-326

Sealing machine



Order code:AM-327

Sealing machine



Order code:AM-328

Sealing machine



Order code:AM-329

Sealing machine



Order code:AM-330

Dental portable printer



Order code:AM-331

heat Sealing sterilization reels



Order code:AM-332

Sterillization pouch



Order code:AM-333

Composite resin heater



Order code:AM-334

Composite resin heater



Order code:AM-335

Composite resin heater



Order code:AM-336

Dental lab equipment

Plastic Water Distiller



Order code:AM-337

Stainless water distiller



Intelligent version

Order code:AM-338

Stainless water distiller



One button version

Order code:AM-339

Stainless water distiller



Two buttons version

Order code:AM-340

Stainless Steel High Pressure Pot



Order code:AM-341

Serilier



12L

Order code:AM-342

Serilier



18L

Order code:AM-343

Serilier



24L

Order code:AM-344

Vacuum molding machine



Order code:AM-345

Dental Vibrator




Order code:AM-346

Magnetic Dental Model Articulator



Order code:AM-347

Dental Lab Digital Wax Pot



Order code:AM-348

Dental Lab Digital Wax Pot



Order code:AM-350

Dental Lab Digital Wax Pot



Order code:AM-351

Amalgamator



Order code:AM-352

Amalgamator




Order code:AM-353

Dental tools



Handpiece repair tool



Camera photography flash I




Dental Oral LED Light



Silicone Dispenser Gun



Composite filling instruments



Aesthetic Plugger



Dental tool set



Plastic teeth Dental
Impression Tray



Stainless teeth Dental
Impression Tray



10pcs Extraction forceps sets



13pcs Dental root elevators set



Dental Orthodontic Crowns
Removal Complete Set



Dental Forming Sheet Clip



dental instrument
Measuring Loop



Oral Planting Positioning
Guide Dental Implant
Surgery



Stainless Steel Dental
Implant Osteotomes
Instruments



Dental Iwanson Gauge



stainless syringe



Ultrasonic scaler tips holder



Orthodontic Instruments Dental Elastic
Separating Pliers



Dental Finger Ruler



Dental Mechanic Screw
Driver



Dental Rubber Dam Clamps



Dental Broken Root Drill
Remnant



Hydraulic Press Polymerization Unit

Disposable material



10pcs high speed dental diamond bur set



30pcs high speed dental diamond bur set



Porcelain teeth polishing kits



Compsite polishing kits



Dental Implants Orthodontic Micro Screw Stainless Steel(120pcs /240pcs)



Dental Implants Orthodontic Micro Screw Gold Plated(120pcs /240pcs)



Dental Implants Orthodontic Micro Screw Titanium(120pcs)



Disc Rubber Polisher



Dental carbide drill



Dental bur box



Russian Polishing Discs




H/K Files



triple way syringe nozzle tips



Colorful Instrument Code Rings



Endo files



dental ceramic mounted stones



Low Speed Handpiece Polisher



Dental Laboratory Ceramics Cutting



Gloves



Disposable applicator



Disposable Propphy Polishing Brush



Diamond Stones Zirconia Polishing(Blue)



Diamond Stones Zirconia Polishing(Green)



Wheel Dental Polishing Disc



Polishing strips



Dental Orthodontic Elastic Power Chains



dental Rubber Dam Clamp Resin Clip



Sectional matrix band

Disposable material



Sectional matrix band



Impression Material



Articulating Paper



Sectional matrix band



Disposable Matrix Band



Colorful Dental Wax
Orthodontic



Rubber pure latex dental
dam



Remove tools



Irrigation needle tips



Microdont Dental Metal
Polishing Stick Strip



Disposable wooden wedges



Disposable Plastic Wedges



Dental Guard Wedges with
metal



Polishing powder



dental needles



Dental suction tips



Dental Syringe



T shape cheek retractor



Impression mixing tips



Cotton Roll Dispenser



Elastic Ligature Tie



High speed handpiece 7balls
steel bearing



Tooth model



Tooth model



Teeth bracket



Niti and stainless arch wires



lingual retainers



Buccal tube



Orthodontic Dental Bracket
Gauge Locator



Dental Orthodontic Niti Open Coil Spring



Traction hook



Open surgical curved long hook



Crimpable hook Long



Crimpable stop